IVORY DENTIN GRAFT™ – INSTRUCTIONS FOR USE – SYRINGE

1 CAUTION ⚠

CAUTION: Federal law restricts this device to sale by or on the order of a dentist or physician.

Carefully read all instructions prior to use.

Failure to observe warnings and precautions noted throughout these instructions may result in complications. Any recommendations within and are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

2 INDICATION

Intended Use:

Ivory Dentin Graft™ is a medical device intended to be used as a bone graft material for the repair or augmentation of bone defects in dental procedures.

Indications for Use with recommended quantity of Ivory Dentin Graft™ to be applied on the bone defect:

Indication	Recommended quantity (g)*
1. Augmentation or reconstructive treatment of alveolar ridge	Varies
2. Filling of infra-bony periodontal defects	0.25-1.0
3. Filling of defects after root resection, apicoectomy, and cystectomy	0.25-1.0
4. Filling of extraction sockets to enhance preservation of the alveolar ridge	0.5–2.0
5. Elevation of maxillary sinus floor	2.0-5.0
6. Filling of periodontal defects in conjunction with products intended for Guided Tissue Regeneration (GTR) and Guided Bone Regeneration (GBR)	0.25–1.5
7. Filling of peri-implant defects in conjunction with products intended for Guided Bone Regeneration (GBR)	0.25–1.5

* Other clinical parameters may also impact the recommended quantity to be administered such as: Size and anatomical site of the defect, Bone quality at the recipient site, Thickness of the bone, Poor blood supply, Maxillary or Mandibular defect, Type of teeth (e.g. incisors, molars), Functional load at defect for a daily dental activity requiring, Pressure from adjacent tissues, External or internal augmentation procedure.

3 DESCRIPTION

Ivory Dentin Graft™ is a porcine sourced, biocompatible bone graft material obtained using standardized, controlled manufacturing processes.

It consists of porous granules of hydroxyapatite which retain the natural form of the source porcine dentin and also the natural protein matrix which consists largely of porcine collagen. Ivory Dentin Graft™ undergoes rigorous cleaning during manufacture and is sterilized by gamma-irradiation.

4 CONTRAINDICATIONS

- Ivory Dentin Graft™ is contraindicated under the following conditions:
 - In the presence of acute inflammation.
 - In patients with known or suspected hypersensitivity, or allergy, to porcine products.
 - In patients affected by immunological disorders.
- The following are relative contraindications:
 - Serious bone diseases of endocrine etiology.
 - Serious disturbances of bone metabolism.
 - Ongoing treatment with glucocorticoids, mineralocorticoids or with agents affecting calcium metabolism (e.g. calcitonin).
 - Severe or difficult to control diabetes mellitus.
 - Irradiation therapy, chemotherapy or immunosuppressive therapy in the last 5 years.
 - Malignancies (because the value of diagnostic X-ray examinations in case of a tumor recurrence at the site of the implant is reduced).

5 WARNINGS & PRECAUTIONS ⚠

- Do not use product that has damaged or opened packaging, as sterility may be compromised, use of non-sterile product may result in patient injury.
- Do not use product that has expired.
- Ivory Dentin Graft™ is supplied sterile for single use only. The device should not be re-sterilized. Re-sterilizing the product may result in patient injury. Do NOT attempt to re-use the device for more than one patient, as this can cause cross-contamination or infection.
- Ivory Dentin Graft™ should be carefully examined prior to surgery and continuously monitored throughout the surgical procedure to ensure the structural integrity and sterility of the device has not been compromised in any way. Using a damaged device may result in patient injury.
- If the product is removed from the patient, it should be discarded in accordance with recognized procedures to discard regulated medical waste materials.
- Effect on pediatric patients is not known.
- Effect on pregnant or lactating women is not known.
- MR Magnetic resonance statement: Ivory Dentin Graft™ has not been evaluated for safety and compatibility in MRI environment.
- Any serious incident that has occurred in relation to the device should be reported immediately to Ivory Graft Ltd. and the competent authority of the Member State in which the user and/or patient is established.

6 ADVERSE REACTIONS

As Ivory Dentin Graft™ is bone graft material, allergic reactions may not be totally excluded.

Possible allergic reaction due to porcine cellular remnants, such as Alpha-Gal syndrome.

Possible complications which may occur with any surgery include swelling at the surgical site, flap sloughing, bleeding, dehiscence, hematoma, increased sensitivity and pain, bone loss, redness, fibrous adhesions, local inflammation, failure of implant osseo-integration and bone formation requiring additional grafting.

7 DIRECTION FOR USE

The general principles of sterile handling and patient medication must be followed when using Ivory Dentin Graft™.

Preparation: Ivory Dentin Graft™ should be handled using sterile gloves and/or sterile instruments. The filling of bony defects with Ivory Dentin Graft™ requires appropriate preparation of the defect. After exposure of the defect, the defect walls should be debrided and all granulation tissue removed. The particles should be placed in direct contact with well vascularized, bleeding bone surfaces. Cortical bone should be mechanically perforated to facilitate ingrowth of new blood vessels and bone forming cells. Ensure that the product can be properly stabilized and protected at the graft site. Placement of a membrane over the particulate bone graft would be according to the physician decision based on each particular case.

- Before injecting Ivory Dentin Graft™, remove the Net from the Syringe tip using sterile forceps. Keep the net in the sterile area in case it needs to be re-applied..
- Immediately apply Ivory Dentin Graft™ in areas where the graft can be adequately contained. The defect should be completely filled ensuring good contact with the walls, however excessive pressure or packing of the material should be avoided as this may inhibit bone growth into the graft. In case of granules jamming, re-apply the net and draw up additional saline or blood into the Syringe and gently tap to loosen granules, expel excess liquid, remove the net, as before and apply to the surgical site.
- Do not compromise blood supply to the defect area.
- Modelling of the graft material using a sterile spatula or other suitable instrument may be performed after application.
- Overfilling of defects is to be avoided.
- Moderate pressure may be applied to reduce bleeding.
- The grafted particles must be packed into a stable form by applying constant gentle pressure. The area is then sutured to close the gum tissue. According to the physician decision a membrane barrier, commonly used in dental practice, may be used to cover the granules in order to prevent granule spread.
- Dispose of any unused product in accordance with recognized procedures to discard regulated medical waste materials.

Site closure:

- a) Wound closure should ensure complete covering of the graft material, using primary flap mobilization if necessary. The flap should be sutured to achieve further closure without any tension. The defect should be fully closed with no gaping remaining.
- b) Double suturing is strongly recommended due to the anticipated postoperative swelling.
- c) If particles are displaced into the surrounding soft tissue during the grafting procedure, these particles should be removed.
- d) Apply a surgical dressing on the surgical area as needed for 1 to 2 weeks.

8 POST-OPERATIVE CARE

- The Dentist should instruct patients not to apply pressure to the surgical site during the bone healing process.
- The Dentist should inspect and evaluate patients for any allergic response due to the graft.
- The graft site should be allowed to heal for at least 3 to 6 months prior to implant placement.
- Successful periodontal treatment requires adequate oral hygiene and infection control. It is recommended to have a hygiene phase, including instruction to the patient, both prior to and following the procedure.
- Granules removal is seldom required. However, in case of severe systemic infection or traumatic injury, early graft removal may be necessary. In this case, the graft site should be surgically opened and all the graft granules together with the adjacent tissue should be thoroughly removed by curettage.

9 HOW SUPPLIED /SPECIFICATIONS

Ivory Dentin Graft™ is supplied sterile in a pouch packaged syringe. Different quantities are available per package in order to match the clinical requirements:

Article Number	Product	Weight	Particle Size
IV-050-S	Ivory Dentin Graft™ Syringe	0.50g	300-900µm
IV-100-S		1.00g	300-900µm

10 STORAGE CONDITIONS

- Store protected from direct sunlight or contact with hot surfaces, in a dry, clean and ventilated place at a temperature between +5°C - +30°C /+1°F- 86°F.

11 EXPLANATION OF SYMBOLS

Symbol	Explanation	Symbol	Explanation
	Qty. (in grs.) per package		Do not re-sterilize
	The number of the NB next to CE marking		Do not use if package is damaged
	Manufacturer		Do not re-use
	Date of manufacture		Keep away from sunlight
	Authorized representative in the European Community		Keep dry
	Importer		Consult instructions for use
	Distributor		Caution
	Use by date		Non-pyrogenic
	Batch code		Contains biological material of animal origin
	Catalogue number		Medical device
	Sterilized using irradiation		Temperature limit
	Federal US law restricts this device for sale by or only on the order of a physician		

12 GENERAL INFORMATION

- Ivory Graft Ltd., 1 Ha Tahana Street, Kefar Sava 4453001, ISRAEL. Website: www.ivorygraft.com
- Email: info@ivorygraft.com
- Emergo Europe B.V., Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands
- MedEnvoy, Prinses Margrietplantsoen 33 - Suite 123, 2955 AM The Hague, The Netherlands
- The SSCP is available in the European Database on Medical Devices (EUDAMED), which can be accessed in the following website: https://ec.europa.eu/tools/eudamed, using the Basic UDI-DI: 7290018472JGD001VG
- Date of Issue: 12/12/2022

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- Ivory Graft Ltd., 1 HaTahana Street, Kefar Sava 4453001, ISRAEL. Website: www.ivorygraft.com
- Emergo Europe B.V., Westervoortsedijk 60, 6827 AT Arnhem, Pays-Bas